

DispensMate Bottle Top Dispenser
(Brown reagent bottle is optional)



- Excellent chemical resistance, components are made of PTFE, FEP, BSG, PP
- Fully autoclavable at 121°C
- Four ranges of bottle-top dispenser covering a volume range from 0.5mL to 50mL
- Easy for cleaning and maintenance
- The optional flexible discharge tube with safety handle permits fast and

precise dispensing

- Vapor pressure max. 500mbar, viscosity max. 500mm²/s,
- temperature max. 40oC, density max. 2.2g/cm³
- DispenseMate is supplied with S40, GL32, GL38, GL25, GL28

Cat.No.	Volume Range ml	Graduation ml	A $\leq$ $\pm$		CV $\leq$	
			%	μ l	%	μ l
7032100001	0.5-5	0.1	0.5	25	0.1	5
7032100002	1.0-10	0.2	0.5	50	0.1	10
7032100003	2.5-25	0.5	0.5	125	0.1	25
7032100004	5.0-50	1.0	0.5	250	0.1	50