

DZ Series C-type stainless steel electric double distilled water device

DZ5C/DZ10C/DZ20C



Performance and characteristic:

Machine is made of high quality stainless steel. Twin towers, reasonable structure. Condensing coil structure, the cooling effect is good, large volume of water. Dry automatic cut off the heating power, water to be added back to arrive at the job level and continue heating. Initial steaming and re-distilled water control were set. Copper immersion electric heating pipe, high efficiency, long service life.

Main Technical Data:

Model	Electric Supply	Power	Water discharge	Product size(Wooden)
DZ5C	380V/50Hz	10KW	≥ 5L/h	42*72*86
DZ10C	380V/50Hz	17KW	≥ 10L/h	48*75*100
DZ20C	380V/50Hz	32KW	≥ 20L/h	56*87*120