

HH-S6 Digital Thermostatic Water Bath



Features

Constant Temperature Water Bath

1. High quality stainless steel material
2. Heat resisting and anticorrosion
3. Operational temperature controller and bi-integral A/D converter
4. Digital display,durable

Specification

Specifications of HH-S6 Digital Constant Temperature Water Bath	
Power Supply	AC 220V/110V 50Hz/60Hz
Temperature control range	Room temperature to 100 Celsius
Temperature fluctuation	0.5 Celsius
Heating power	1500W
Working room size	L50 W30H15cm
Digital display	Yes
Laboratory water bath of easy to use, quality stainless steel material, good reliability.	

Option

1. With one whole cover
- 2.High temperature controller with accuracy $\pm 0.01^{\circ}\text{C}$

