

MicroPette Plus



- Fully autoclavable
- Ergonomic design provides excellent operating experience
- Easy-to-read volume display
- The pipettes cover the volume range from 0.1µL to 5mL
- Easy calibration and maintenance
- Manufactured from innovative materials
- Each MicroPette Plus supplied with an
- individual calibration certificate according to ISO8655

Volume Range	Increment	Test Volume (µl)	Accuracy error		Precision error	
			%	µl	%	µl
0.1-2.5µl	0.05µl	2.5	2.50%	0.0625	2.00%	0.05
		1.25	3.00%	0.0375	3.00%	0.0375
		0.25	12.00%	0.03	6.00%	0.015
0.5-10µl	0.1µl	10	1.00%	0.1	0.80%	0.08
		5	1.50%	0.075	1.50%	0.075
		1	2.50%	0.025	1.50%	0.015

2-20µl	0.5µl	20	0.90%	0.18	0.40%	0.08
		10	1.20%	0.12	1.00%	0.1
		2	3.00%	0.06	2.00%	0.04
5-50µl	0.5µl	50	0.60%	0.3	0.30%	0.15
		25	0.90%	0.225	0.60%	0.15
		5	2.00%	0.1	2.00%	0.1
10-100µl	1µl	100	0.80%	0.8	0.15%	0.15
		50	1.00%	0.5	0.40%	0.2
		10	3.00%	0.3	1.50%	0.15
20-200µl	1µl	200	0.60%	1.2	0.15%	0.3
		100	0.80%	0.8	0.30%	0.3
		20	3.00%	0.6	1.00%	0.2
50-200µl	1µl	200	0.60%	1.2	0.15%	0.3
		100	0.80%	0.8	0.30%	0.3
		50	1.00%	0.5	0.40%	0.2
100-1000µl	5µl	1000	0.60%	6	0.20%	2
		500	0.70%	3.5	0.25%	1.25
		100	2.00%	2	0.70%	0.7
200-1000µl	5µl	1000	0.60%	6	0.20%	2
		500	0.70%	3.5	0.25%	1.25
		200	0.90%	1.8	0.30%	0.6
1000-5000µl	50µl	5000	0.50%	25	0.15%	7.5
		2500	0.60%	15	0.30%	7.5
		1000	0.70%	7	0.30%	3
2-10ml	0.1ml	10ml	0.60%	60	0.20%	20
		5ml	1.20%	60	0.30%	15
		2ml	3.00%	60	0.60%	12