

SP-100/200/300 Multi-Positon Magnetic Stirrer

SP-100 / SP-200 / SP-300 series multi-position magnetic stirrer is designed for the high-efficiency experiments and has 4/8/12 stir bar configurations to choose. It can use laboratory space effective. DC brushless motor make the model:low speed but stable,high speed but powerful with low noise. Stainless steel is coated by silicone resin film disk. Non-slip and anti-corrosion.



Speed range adjustable



Strong Magnetic



DC Brushless Motor



Time Setting



【 Product Features 】

1. Digital display stirring speed at real-time reaction conditions.
2. Efficient multi-position magnetic stirrer configure 4/8/12 stirring positions to use lab space effective.
3. DC brushless motor.Low speed but stable,high speed but powerful.Low noise and free maintenance.
4. Stainless steel is coated by silicone resin film disk. Non-slip and anti-corrosion.

【 Technical Parameters 】

Model	SP-100	SP-200	SP-300
Stirring Point Quantity	4	8	12
Max. Size of Stirrer Stick	40mm	40mm	40mm
Speed Range	200 ~1200rpm	200 ~1200rpm	200 ~1200rpm
Platform Material	stainless steel with silicon pad	stainless steel with silicon pad	stainless steel with silicon pad
Motor	brushless DC motor	brushless DC motor	brushless DC motor
Voltage	AC100 ~230V, 50/60Hz	AC100 ~230V, 50/60Hz	AC100 ~230V, 50/60Hz
Power	20W	20W	20W
Fuse	250V,1A, $\Phi 5 \times 20$	250V,1A, $\Phi 5 \times 20$	250V,1A, $\Phi 5 \times 20$
Dimension (WxDxH)	115 x 480 x 50mm	205 x 480 x 50mm	295 x 480 x 50mm
Net Weight	3kgs	4kgs	5kgs